
INSTRUMENT IDENTIFICATION

INNOVATIVE LAB SERVICES LLC - API 4000 LC/MS/MS SYSTEM - PLANNED MAINTENANCE
PROCEDURE (Doc 1_1, Rev 1.1)

Instrument Serial #

Eyeye

Instrument Name / INV #

Test

CUSTOMER INFORMATION

Customer Site	Nestle
Serial Number	2646464
Customer Name	Nextle
Service Engineer	Harris
Date	07/23/2026

SIGN OFF

By signing this form you acknowledge that the planned maintenance was completely performed.
Customer _____

Service Engineer _____



PRE-PM 1.0 VACUUM PRESSURE CHECK

Guideline: CAD = 0 -> $\leq 0.8 \times 10^{-5}$ torr (0.08 V). CAD = 12 -> $\leq 5.5 \times 10^{-5}$ torr (+/- 0.3×10^{-5}).

Vacuum pressure, CAD = 0 _____ .25

Vacuum pressure, CAD = 12 _____

PRE-PM 2.0 POSITIVE ION MODE TEST (10 SCANS MCA)

Guideline (10 Scans MCA): Q1 m/z 906.673 $\geq 2.0 \times 10^7$ | Q1 m/z 2242.6 $\geq 1.0 \times 10^6$ | Q3 m/z 906.673 $\geq 2.0 \times 10^7$ | Q3 m/z 2242.6 $\geq 8.0 \times 10^5$. FWHM 0.6 - 0.8 amu throughout.

Q1 m/z 906.673

Intensity _____

FWHM _____

Q1 m/z 2242.6

Intensity _____

FWHM _____

Q3 m/z 906.673

Intensity _____

FWHM _____

Q3 m/z 2242.6

Intensity _____

FWHM _____

PRE-PM 3.0 NEGATIVE ION MODE TEST (10 SCANS MCA)

Guideline (10 Scans MCA): Q1 m/z 933.636 $\geq 1.1 \times 10^7$ | Q1 m/z 2037.4 $\geq 1.0 \times 10^6$ | Q3 m/z 933.636 $\geq 9.0 \times 10^6$. FWHM 0.6 - 0.8 amu throughout. Q3 m/z 2037.4 is N/A.

Q1 m/z 933.636

Intensity

FWHM

Q1 m/z 2037.4

Intensity

FWHM

Q3 m/z 933.636

Intensity

FWHM

PRE-PM 4.0 RESERPINE MS/MS TEST (10 SCANS MCA)

Guideline (10 Scans MCA): MS/MS transmission efficiency (m/z 195.0 / attenuated m/z 609.3) x 100 ≥ 10.0 %

Attenuated m/z 609.3

Intensity

FWHM

0.70

m/z 195.0

Intensity

FWHM

Transmission Efficiency %

PM SECTION - MAINTENANCE AND CLEANING

1. Check operation of the Turbo V ion source heaters ☒
2. Check lens voltages (Lens Power Supply module) ☒
3. Check RF tuning voltages (Q1 and Q3 QPS AMP module). Only if required, tune QPS AMPs. ☒
4. Check operation of the instrument cooling fans ☒
- CEM voltage (V)
- SHUTDOWN AND PERFORM MAINTENANCE / CLEANING
- AC input voltage (V)
7. Inspect gas tubing for breakage (CAD gas, interface gas tubes, etc) ☐
8. Remove and clean the curtain plate ☐
9. If necessary, remove and clean the orifice plate ☐
10. If necessary, clean the front end of Q0 ☐
11. Replace the air filter ☐
12. Replace the roughing pumps oil ☐
13. Inspect the oil return systems and if necessary replace the oil filter ☐
14. If applicable, perform PM on Peak Scientific gas generator ☐
15. Reinstall all covers and pump down ☒
16. Re-check operation of the instrument cooling fans ☒
17. If necessary, perform Turbo V ion source maintenance ☒
18. Clear computer temp files (Start / Search / Files and Folders -> *.tmp -> Local Hard Drive (C), (D) -> Search now -> Edit / Select All / Files / Delete / Empty Recycle Bin) ☐

PERFORMANCE VERIFICATION 1.0 GAS PRESSURE

Low Gas 1 / Gas 2 pressure can affect turbo spray and sensitivity. Source Exhaust pressure above 60 psi can cause exhaust valve failures.

Requirement: Gas 1 / Gas 2 100 - 105 psi | Curtain / CAD Gas 60 psi | Source Exhaust Supply 55 - 60 psi

Gas 1 / Gas 2 (psi)

Curtain / CAD Gas (psi) 60.0

Source Exhaust Supply (psi)

PERFORMANCE VERIFICATION 2.0 VACUUM PRESSURE

Guideline: CAD = 0 -> $\leq 0.8 \times 10^{-5}$ torr (0.08 V). CAD = 12 -> $\leq 5.5 \times 10^{-5}$ torr (+/- 0.3×10^{-5}).

Vacuum pressure, CAD = 0

Vacuum pressure, CAD = 12

PERFORMANCE VERIFICATION 3.0 POSITIVE ION MODE TEST (10 SCANS MCA)

Guideline (10 Scans MCA): Q1 m/z 906.673 $\geq 2.0 \times 10^7$ | Q1 m/z 2242.6 $\geq 1.0 \times 10^6$ | Q3 m/z 906.673 $\geq 2.0 \times 10^7$ | Q3 m/z 2242.6 $\geq 8.0 \times 10^5$. FWHM 0.6 - 0.8 amu throughout.

Q1 m/z 906.673

Intensity

FWHM 0.60

Q1 m/z 2242.6

Intensity

FWHM

Q3 m/z 906.673

Intensity

FWHM

Q3 m/z 2242.6

Intensity

FWHM

PERFORMANCE VERIFICATION 4.0 NEGATIVE ION MODE TEST (10 SCANS MCA)

Guideline (10 Scans MCA): Q1 m/z 933.636 $\geq 1.1 \times 10^7$ | Q1 m/z 2037.4 $\geq 1.0 \times 10^6$ | Q3 m/z 933.636 $\geq 9.0 \times 10^6$. FWHM 0.6 - 0.8 amu throughout. Q3 m/z 2037.4 is N/A.

Q1 m/z 933.636

Intensity

FWHM

Q1 m/z 2037.4

Intensity

FWHM

Q3 m/z 933.636

Intensity

FWHM

PERFORMANCE VERIFICATION 5.0 RESERPINE MS/MS TEST (10 SCANS MCA)

Guideline (10 Scans MCA): MS/MS transmission efficiency (m/z 195.0 / attenuated m/z 609.3) x 100
 $\geq 10.0 \%$

Attenuated m/z 609.3

Intensity

FWHM

m/z 195.0

Intensity

FWHM

Transmission Efficiency %

6.0 TEST SPECTRA AND COMMENTS

Test spectra printed and attached

☐

Comments or Observations

This is a test

Annotation

This is a test